

Bis-(Mal-PEG3)-PH-N-succinimidyl acetate

Chemical Properties

CAS No. :

Formula: C₃₁H₃₅N₃O₁₄

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bis-(Mal-PEG3)-PH-N-succinimidyl acetate is a PEG-based PROTAC linker used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481