

Mal-PEG12-NH-Boc

Chemical Properties

CAS No. : 2170654-80-3

Formula: C₃₅H₆₄N₂O₁₆

Molecular Weight: 768.89



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Mal-PEG12-NH-Boc is a PROTAC linker classified under the PEG category and is used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3006 mL	6.5029 mL	13.0058 mL
5 mM	0.2601 mL	1.3006 mL	2.6012 mL
10 mM	0.1301 mL	0.6503 mL	1.3006 mL
50 mM	0.026 mL	0.1301 mL	0.2601 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481