

Fmoc-PEG9-NHS ester

Chemical Properties

CAS No. :

Formula: C₄₀H₅₆N₂O₁₅

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fmoc-PEG9-NHS ester is a type of PROTAC linker that belongs to the PEG category. It is used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481