

## Propargyl-PEG11-acid

## Chemical Properties

CAS No. :

Formula: C<sub>26</sub>H<sub>48</sub>O<sub>13</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Propargyl-PEG11-acid is a PROTAC linker belonging to the PEG class, suitable for synthesizing PROTAC molecules. As a click chemistry reagent, it contains an alkyne group that can undergo a copper-catalyzed azide-alkyne cycloaddition reaction (CuAAC) with molecules containing azide groups.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481