

Thalidomide-NH-amido-C3-NH2

Chemical Properties

CAS No. : 2763827-11-6


Formula: C18H21N5O5

Molecular Weight: 387.39

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Thalidomide-NH-amido-C3-NH2 is a synthetic E3 ligase ligand-linker conjugate, consisting of a Thalidomide-based cereblon ligand and a linker, designed for the synthesis of PROTAC.
Targets(IC50)	Apoptosis, Autophagy, E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5814 mL	12.9069 mL	25.8138 mL
5 mM	0.5163 mL	2.5814 mL	5.1628 mL
10 mM	0.2581 mL	1.2907 mL	2.5814 mL
50 mM	0.0516 mL	0.2581 mL	0.5163 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481