

Thalidomide-NH-amido-C3-NH2 hydrochloride

Chemical Properties

CAS No. :

Formula: C18H22ClN5O5

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Thalidomide-NH-amido-C3-NH2 hydrochloride is a synthetic E3 ligase ligand-linker conjugate, incorporating a Thalidomide-based cereblon ligand and a linker, which is useful for PROTAC synthesis. |
| Targets(IC50) | Apoptosis, Autophagy, E3 Ligase Ligand-Linker Conjugates |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481