

## Pomalidomide-5-C2-NH2 hydrochloride

### Chemical Properties

CAS No. :

Formula: C<sub>15</sub>H<sub>17</sub>ClN<sub>4</sub>O<sub>4</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Pomalidomide-5-C2-NH2 (hydrochloride) is an E3 ligase ligand-linker conjugate used in synthesizing PROTACs.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481