

## MS83 epimer 1

## Chemical Properties

CAS No. :

Formula: C<sub>52</sub>H<sub>55</sub>ClN<sub>10</sub>O<sub>8</sub>S<sub>2</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	MS83epimer 1 is the epimer of MS83, a proof-of-concept PROTAC compound. MS83 is created by linking a KEAP1 ligand to a BRD4/3/2 binder.
Targets(IC50)	PROTACs

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481