

## Aminoxy-PEG3-bromide hydrochloride

### Chemical Properties

CAS No. :

Formula: C<sub>8</sub>H<sub>19</sub>BrClNO<sub>4</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Aminoxy-PEG3-bromide hydrochloride is a PEG-based PROTAC linker used in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481