

TCO-PEG2-Sulfo-NHS ester sodium

Chemical Properties

CAS No. :

Formula: C₂₀H₂₉N₂NaO₁₁S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TCO-PEG2-Sulfo-NHS ester sodium is a PROTAC linker from the PEG category, utilized in the synthesis of PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481