

Boc-L-Cys(Propargyl)-OH (DCHA)

Chemical Properties

CAS No. :

Formula: C₂₃H₄₀N₂O₄S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Boc-L-Cys(Propargyl)-OH (DCHA) is a reagent utilized in click chemistry and can be applied in various biochemical studies.
Targets(IC50)	ADC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481