

Propargyl-PEG-amine

Chemical Properties

CAS No. : 956348-61-1
Formula: $(C_2H_4O)_n C_5H_9NO$
Molecular Weight:
Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Propargyl-PEG-amine is a PROTAC linker within the PEG class, suitable for synthesizing PROTAC molecules.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481