

dTAG-47-NEG

Chemical Properties

CAS No. :

Formula: C59H75N5O13

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	dTAG-47-NEG is an analogue of dTAG-47 that is unable to bind and recruit cereblon (CRBN), serving as a heterobifunctional negative control for dTAG-47.
Targets(IC50)	PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481