

Pomalidomide-C3-NHS ester

Chemical Properties

CAS No. : 2828476-87-3


Formula: C₂₁H₂₀N₄O₈

Molecular Weight: 456.41

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pomalidomide-C3-NHS ester is a synthetic E3 ligase ligand-linker conjugate that comprises a Pomalidomide-based cereblon ligand coupled with an NHS ester.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.191 mL	10.9551 mL	21.9101 mL
5 mM	0.4382 mL	2.191 mL	4.382 mL
10 mM	0.2191 mL	1.0955 mL	2.191 mL
50 mM	0.0438 mL	0.2191 mL	0.4382 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481