

Cathepsin Inhibitor 4

Chemical Properties

CAS No. :

Formula: C₂₄H₃₅N₃O₅

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cathepsin Inhibitor 4 (Compound 45) is a selective inhibitor of Cathepsin S, with a Ki value of 0.04 nM.
Targets(IC50)	Cysteine Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481