


(S)-OY-101

Chemical Properties

CAS No. : 2924484-79-5
 Formula: C₂₇H₃₁N₂O₄
 Molecular Weight: 433.54
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-OY-101 serves as the isomer of OY-101.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3066 mL	11.533 mL	23.0659 mL
5 mM	0.4613 mL	2.3066 mL	4.6132 mL
10 mM	0.2307 mL	1.1533 mL	2.3066 mL
50 mM	0.0461 mL	0.2307 mL	0.4613 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481