

## Thalidomide-4-O-C9-NH2 hydrochloride

### Chemical Properties

CAS No. :

Formula: C<sub>22</sub>H<sub>30</sub>ClN<sub>3</sub>O<sub>5</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Thalidomide-4-O-C9-NH2 hydrochloride is a Thalidomide-based cereblon ligand capable of recruiting CRBN protein. This compound can be linked to a target protein ligand via a linker to form PROTAC molecules, such as [THAL-SNS-032].
Targets(IC50)	Ligands for E3 Ligase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481