


PROTAC SMARCA2 degrader-1

Chemical Properties

CAS No. :	2892523-73-6
Formula:	C49H62N10O6S
Molecular Weight:	919.15
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTACSMARCA2 degrader-1 (compound ex7) is a potent BRM (SMARCA2) degrader with a DC50 of less than 0.1 μ M. It functions as a linker for E3 ubiquitin ligase, making it useful in cancer research, and is a conjugate of a target protein ligand and linker.
Targets(IC50)	Target Protein Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.088 mL	5.4398 mL	10.8796 mL
5 mM	0.2176 mL	1.088 mL	2.1759 mL
10 mM	0.1088 mL	0.544 mL	1.088 mL
50 mM	0.0218 mL	0.1088 mL	0.2176 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481