

(S)-Phe-A110/B319

Chemical Properties

CAS No. :

Formula: C62H57F3N8O14S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(S)-Phe-A110/B319 is a selective binder with a 20-fold higher affinity for the H1047R mutant of the p110 α subunit within the p110 α /p85 α PI3K complex.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481