

## Antifungal agent 88

### Chemical Properties

CAS No. :

Formula: C<sub>25</sub>H<sub>22</sub>N<sub>4</sub>O<sub>5</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Antifungalagent 88 exhibits antifungal activity, with a MIC value ranging from <0.125 µg/mL to 1 µg/mL.
Targets(IC50)	Antibacterial

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481