

P-gp inhibitor 16

Chemical Properties

CAS No. :

Formula: C₃₅H₃₅N₅O₄

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	P-gp inhibitor 16 (compound 14) is an inhibitor of p-glycoprotein. It significantly enhances apoptosis induced by doxorubicin and displays anticancer properties.
Targets(IC50)	Apoptosis,P-gp

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481