

MDM2-IN-23

## Chemical Properties

CAS No. :

Formula: C<sub>26</sub>H<sub>20</sub>Cl<sub>4</sub>N<sub>2</sub>O<sub>2</sub>S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	MDM2-IN-23 is an inhibitor of MDM2, exhibiting an IC <sub>50</sub> of 60.09μM against MCF-7 cells.
Targets(IC <sub>50</sub> )	Ligands for E3 Ligase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481