

MRS7935

Chemical Properties

CAS No. :

Formula: C₂₄H₂₀ClF₂NOS

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MRS7935 (compound 15) is a positive allosteric modulator (PAM) of the A1 adenosine receptor (A1AR), with an EC50 of approximately 2 µM.
Targets(IC50)	Adenosine Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481