

## Biotin-PEG3-CoenzymeA

### Chemical Properties

CAS No. :

Formula: C46H84N15O24P3S2

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Biotin-PEG3-CoenzymeA is a non-cell-permeable substrate based on biotin with a PEG linker. It can be utilized for biotinylation in living cells, labeling in solutions, and interaction studies.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481