

(R)-TCO4-PEG2-NH2

Chemical Properties

CAS No. :

Formula: C15H28N2O4

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(R)-TCO4-PEG2-NH2 is a PROTAC linker categorized under the PEG class. It contains a TCO group capable of undergoing a specific "click" reaction with tetrazine groups.
Targets(IC50)	PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481