

## Pomalidomide-amino-PEG4-C4-Cl

## Chemical Properties

CAS No. : 2924835-36-7


Formula: C<sub>27</sub>H<sub>36</sub>ClN<sub>3</sub>O<sub>9</sub>

Molecular Weight: 582.04

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Pomalidomide-amino-PEG4-C4-Cl (Intermediate 4d) is an E3 ligase ligand-linker conjugate featuring a Pomalidomide-based cereblon (CRBN) ligand and a linker. This compound is utilized in the synthesis of PROTACs.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7181 mL	8.5905 mL	17.1809 mL
5 mM	0.3436 mL	1.7181 mL	3.4362 mL
10 mM	0.1718 mL	0.859 mL	1.7181 mL
50 mM	0.0344 mL	0.1718 mL	0.3436 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481