

PLGA-PEG-PLGA (1500-1500-1500) (LA/GA 15:1)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PLGA-PEG-PLGA (1500-1500-1500) (LA/GA 15:1) is a matrix material for cancer drug delivery, characterized by a lactic acid to glycolic acid ratio of 15:1. This compound enhances drug bioavailability, efficacy, solubility, drug encapsulation efficiency, and sustained release while minimizing adverse toxicity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481