

NP-5497-KA

## Chemical Properties

CAS No. :

Formula: C<sub>26</sub>H<sub>39</sub>Cl<sub>3</sub>N<sub>4</sub>O<sub>2</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NP-5497-KA is a $\kappa$ -opioid receptor agonist with oral bioavailability.
Targets(IC50)	Opioid Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481