


EN884

Chemical Properties

CAS No. : 2189497-60-5
 Formula: C₁₄H₁₈N₂O
 Molecular Weight: 230.31
 Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	EN884 is a BRD4 degrader that functions through an SKP1 and proteasome-dependent degradation pathway. It is utilized in the synthesis of proteolysis-targeting chimeras (PROTAC).
Targets(IC50)	Epigenetic Reader Domain,Target Protein Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.342 mL	21.7099 mL	43.4197 mL
5 mM	0.8684 mL	4.342 mL	8.6839 mL
10 mM	0.4342 mL	2.171 mL	4.342 mL
50 mM	0.0868 mL	0.4342 mL	0.8684 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481