

Cholic acid-cysteine-cyanuric chloride complex

Chemical Properties

CAS No. :

Formula: C₃₀H₄₄Cl₂N₄O₆S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cholic acid-cysteine-cyanuric chloride complex is a hapten linker molecule, consisting of an antigen, cholic acid, and the covalently bonded cyanuric chloride reactive group.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481