

ZEp

## Chemical Properties

CAS No. :

Formula: C<sub>23</sub>H<sub>30</sub>O<sub>7</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ZEp is an immunogenic hapten derived from zearalenone (ZEN), utilized in immunoassays after the activation of its carboxyl group into the corresponding N-hydroxysuccinimide ester.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481