

ORI-TRN-002

## Chemical Properties

CAS No. :

Formula: C<sub>10</sub>H<sub>7</sub>N<sub>5</sub>O

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ORI-TRN-002 is an AQP4 inhibitor with an IC <sub>50</sub> value of 3 μM. It is applicable in the field of brain edema research.
Targets(IC <sub>50</sub> )	Aquaporin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481