

## AC mAb-X purification resin

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	AC mAb-X purification resin is a specialized Protein A resin. It features a highly rigid agarose matrix and particle sizes ranging from 40 to 120 µm. The ligand is a recombinant alkaline-resistant Protein A (~24KD), and it serves as a purification resin for ADC.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481