

(2S,3R)-Dap-NE hydrochloride

## Chemical Properties

CAS No. :

Formula: C<sub>18</sub>H<sub>29</sub>ClN<sub>2</sub>O<sub>3</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |   |
|---------------|---|
| Description   | (2S,3R)-Dap-NE hydrochloride is an amino acid residue of the pentapeptide Dolastatin 10. Dolastatin 10 inhibits tubulin polymerization and mitosis, exhibiting anticancer properties. |
| Targets(IC50) | Microtubule Associated  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481