

## Gadoteridol pentahydrate

## Chemical Properties

CAS No. :

Formula: C17H38GdN4O12

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Gadoteridol pentahydrate is a gadolinium-based contrast agent utilized for imaging the central nervous system (CNS) in magnetic resonance imaging (MRI).
Targets(IC50)	Others
In vivo	An intracerebral injection of Gadoteridol (1 mM; administered once) improves the spread of rAAV1 within the rat hippocampus.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481