

## Cyanine5 Boc-hydrazide chloride

### Chemical Properties

CAS No. :

Formula: C<sub>37</sub>H<sub>49</sub>ClN<sub>4</sub>O<sub>3</sub>

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cyanine5 Boc-hydrazide (chloride) is a hydrazine-based fluorescent dye that is labeled with Cyanine5 and protected by a Boc group.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481