

Cyanine5 Boc-hydrazide chloride

Chemical Properties

CAS No. :

Formula: C₃₇H₄₉ClN₄O₃

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cyanine5 Boc-hydrazide (chloride) is a hydrazine-based fluorescent dye that is labeled with Cyanine5 and protected by a Boc group.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481