

BODIPY FL Cholesterol

Chemical Properties

CAS No. :

Formula: C₃₆H₅₀BF₂N₃O₂

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BODIPY FLCholesterol is a green fluorescent lipophilic probe. It can be used to study lipid transport processes and lipid-protein interactions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481