

Biotin-11-UTP

Chemical Properties

CAS No. :

Formula: C₂₈H₄₀Li₃N₆O₁₈P₃S

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Biotin-11-UTP is a uridine triphosphate (UTP) labeled with biotin and features a linker consisting of 11 atoms between the biotin and UTP. This compound is widely used for the in vitro biotinylation of transcribed RNA, serving effectively as a labeled RNA probe.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481