

Fluorescein-12-dATP

Chemical Properties

CAS No. :

Formula: C41H41N6O19P3

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fluorescein-12-dATP is a fluorescent dye utilized as a detection probe. It can be employed to label thrombin-binding aptamers through the catalytic polymerization by TdT (terminal deoxynucleotidyl transferase).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481