

Fura-4F AM

Chemical Properties

CAS No. :

Formula: C43H44FN3O24

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fura-4F AM is a fluorescent calcium indicator that permeates cells. Once inside, this probe is hydrolyzed by cellular esterases and becomes an active chelator.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481