

Fluo-4FF AM

## Chemical Properties

CAS No. :

Formula: C<sub>50</sub>H<sub>46</sub>F<sub>4</sub>N<sub>2</sub>O<sub>23</sub>

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Fluo-4FF AM is a cell-permeable fluorescent calcium indicator. As an analog of Fluo-4, Fluo-4FF has lower affinity for calcium, making it suitable for investigating conditions with elevated intracellular calcium levels.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481