

Vari Fluor 350 TSA (200x)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Vari Fluor 350 TSA (200x) (VF 350 Tyramide) is a fluorescent probe within the Vari Fluor TSA series. These fluorescent probes enable the detection method of high-density in situ labeling of target antigens through horseradish peroxidase (HRP). They are widely used in ELISA, IF, and IHC techniques, providing high detection sensitivity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481