

Sulfo-Cy3-NHS disodium

Chemical Properties

CAS No. :

Formula: C₃₆H₄₁N₃Na₂O₁₃S₃

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cy3-NHS disodium is a derivative of Cy3 (Cyanine3) dye that includes sulfonate ions, providing enhanced water solubility. It is capable of labeling biomolecules with amine groups, such as proteins, antibodies, or peptides.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481