

## ATTO 425 streptavidin

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

|               |   |
|---------------|---|
| Description   | ATTO 425 Streptavidin is a derivative of ATTO 425 with streptavidin, exhibiting peak excitation and emission wavelengths of 439/489 nm. |
| Targets(IC50) | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481