

ATTO 594 NHS ester

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ATTO 594 NHS ester is an innovative fluorescent label based on a Rhodamine structure. It features strong absorption, high fluorescence quantum yield, and excellent thermal and photochemical stability, making it suitable for single-molecule detection and high-resolution microscopy. ATTO 594 NHS ester is a derivative of ATTO 594, designed for labeling proteins or antibodies.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481