

GA-PEG5-bromide

Chemical Properties

CAS No. :

Formula: C₄₀H₇₀BrN₃O₁₁

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | GA-PEG5-bromide is a probe used for molecular labeling. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481