

PbS/CdS QDs

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	PbS/CdS QDs (PbS/CdSe QD) serve as fluorescent probes suitable for in vivo imaging within the second near-infrared window.
Targets(IC50)	Others
In vivo	Utilizing the 2nd-NIR fluorescence imaging around 1300 nm with PbS/CdS QDs provides an optimal signal-to-background ratio and minimal blurring for in vivo imaging.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481