

Ag<sub>2</sub>Te QDs

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:   
Keep away from direct sunlight  
Store at -20°C  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Ag <sub>2</sub> Te QDs (Ag <sub>2</sub> Te) serve as an efficient bioprobe for the second near-infrared window (NIR-II), enabling biological imaging with high tissue penetration depth and exceptional spatiotemporal resolution.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481