

Pb:Ag₂Se QDs

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pb:Ag ₂ Se QDs (Pb:Ag ₂ Se) act as an efficient bioprobe for the second near-infrared window (NIR-II), suitable for bioimaging applications that require deep tissue penetration and high spatiotemporal resolution.
Targets(IC50)	Others
In vitro	Pb:Ag ₂ Se QDs exhibit bright NIR-II photoluminescence, alongside high levels of photostability and colloidal stability, as well as commendable biocompatibility.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481